


Instructions for Use
AIRSUPRA™ (ayr soo' prah)
(albuterol and budesonide)
inhalation aerosol, for oral inhalation use
120-inhalation canister

Read this Instructions for Use before you start using AIRSUPRA and each time you get a new refill. There may be new information. This information does not take the place of talking to your healthcare provider about your medical condition or treatment.

Important information

- **For oral inhalation use only.**
- Use AIRSUPRA exactly as your healthcare provider tells you to.
- Before you use AIRSUPRA for the first time, you will need to prime it.
- Rinse the white actuator 1-time each week.
- 2 puffs of medicine is 1-dose. **Do not** use more than 6 doses (12 puffs) in a 24-hour period.
- Rinse your mouth with water, if available, after using AIRSUPRA to decrease the chance of getting a fungal infection (thrush) in your mouth and throat.

Parts of your AIRSUPRA inhaler



Reading the puff indicator

- The puff indicator will count down by 1 each time you spray a puff of medicine.

Pointer


Points to number of puffs remaining.

Yellow zone

Order a new inhaler when the pointer is in the yellow zone.

Red zone

Throw away your inhaler when the pointer is at 0 in the red zone.




Do not try to take a puff of medicine when the pointer is at 0 because you will not receive a full dose.

Before using your AIRSUPRA inhaler

1. Write the Use by date


Write the Use By Date (12 months after you open the foil pouch) on the actuator label using permanent ink.



Write the Use By Date
on the actuator label.


2. Prime your AIRSUPRA inhaler before first use

Take the cover off the mouthpiece by squeezing both sides of the cover while pulling it down.



Shake the inhaler well. Hold your AIRSUPRA inhaler facing away from you and press down firmly on the top of the puff indicator to spray 1 test-puff into the air. Repeat the priming steps for a total of **4 test-puffs**. Shake your AIRSUPRA inhaler before each test-puff.

**Shake before
each test-puff.
Repeat 4 times.**



Extra puffs are provided for priming. **Do not skip priming.**

When and how to re-prime your AIRSUPRA inhaler


Re-prime your inhaler:

- if it is not used for more than 7 days.
- if it is dropped
- after rinsing white actuator

To re-prime your inhaler:

Take the cover off the mouthpiece and spray **2 test-puffs into the air**. Shake your AIRSUPRA inhaler before each test-puff.

**Shake before each
test-puff.
Repeat 2 times.**




Extra puffs are provided for priming. **Do not skip priming.**

How to use your inhaler – Inhaling your medicine

1. Check


Take the cover off the mouthpiece. **Check the mouthpiece for foreign objects and remove objects before use.**




Check inside mouthpiece.

2. Shake and Inhale

Shake well




Breathe out fully



Place the mouthpiece into your mouth and close your lips around the mouthpiece.



Start to breathe in deeply and slowly while spraying **1 puff**. Continue breathing in until you cannot anymore.




Remove the mouthpiece from your mouth. Hold your breath for as long as you can or up to 10 seconds.




1. Remove canister.

Remove the canister and set it aside. Take the cover off the mouthpiece.




2. Rinse through both ends

Rinse each end of the actuator with warm water for 30 seconds. Shake off as much water as you can. **Do not use soap. Do not dry with a towel or tissue.**



60 seconds total



Look into the mouthpiece and actuator to make sure any medicine build-up has been completely washed away. If there is any build-up, repeat the weekly rinse step 2 again.



Check inside mouthpiece and actuator.

3. Repeat

Immediately repeat the shake and inhale step to take a second puff.



2 puffs equals 1 dose.


4. Place the mouthpiece cover back on

Place the mouthpiece cover back on the mouthpiece.




5. Rinse

Rinse your mouth with water, if available (**Do not swallow the water**).



3. Air-dry

Let the actuator air-dry, such as overnight.



Weekly rinse – Rinsing your actuator 1 time each week


- Rinse the white actuator 1-time each week so that the medicine does not build up and block the spray through the mouthpiece.
- **Do not allow the canister to get wet.**
- Re-prime your AIRSUPRA inhaler after rinsing.

4. Re-assemble and re-prime

Re-assemble AIRSUPRA when the actuator is fully dry.

First, place the mouthpiece cover back on. After the mouthpiece cover is on, gently press the canister down into actuator. **Do not insert the canister while the mouthpiece cover is off.**

Re-prime the inhaler by following the Re-prime steps above.



Throwing away your AIRSUPRA inhaler

Throw away your AIRSUPRA inhaler in the household trash:

- when the dose indicator reaches zero “0” or 12 months after you open the foil pouch, whichever comes first.

Do not reuse or use the actuator with medicine canisters from other inhalers.

Do not puncture or throw the canister into a fire or incinerator.

Ordering a new AIRSUPRA inhaler

- Order a new AIRSUPRA inhaler when the pointer on the dose indicator is in the yellow zone.

This Instructions for Use has been approved by the U.S. Food and Drug Administration.

Issued: January 2023

01/23 US-70802 1/23

Emergency inhaler use when wet

If you need to use the inhaler before it is dry:

- shake off as much excess water as you can
- insert the canister into the actuator
- shake and spray 2-test puffs into the air away from your face
- take the prescribed dose
- rinse the actuator again following the Weekly Rinse steps.

Storing your inhaler

- Store your inhaler at room temperature between 68°F to 77°F (20°C to 25°C).
- **Do not store in a humid environment, such as a bathroom.**
- Keep your inhaler and all medicines out of the reach of children.


Instructions for Use
AIRSUPRA™ (ayr soo' prah)
(albuterol and budesonide)
inhalation aerosol, for oral inhalation use
28-inhalation canister

Read this Instructions for Use before you start using AIRSUPRA and each time you get a new refill. There may be new information. This information does not take the place of talking to your healthcare provider about your medical condition or treatment.

Important information

- **For oral inhalation use only.**
- Use AIRSUPRA exactly as your healthcare provider tells you to.
- Before you use AIRSUPRA for the first time, you will need to prime it.
- Rinse the white actuator 1-time each week.
- 2 puffs of medicine is 1-dose. **Do not** use more than 6 doses (12 puffs) in a 24-hour period.
- Rinse your mouth with water, if available, after using AIRSUPRA to decrease the chance of getting a fungal infection (thrush) in your mouth and throat.

Parts of your AIRSUPRA inhaler



Reading the puff indicator

- The puff indicator will count down by 1 each time you spray a puff of medicine.

Pointer


Points to number of puffs remaining.

Yellow zone

Order a new inhaler when the pointer is in the yellow zone.

Red zone

Throw away your inhaler when the pointer is at 0 in the red zone.




Do not try to take a puff of medicine when the pointer is at 0 because you will not receive a full dose.

Before using your AIRSUPRA inhaler


1. Write the Use by date

Write the Use By Date (6 months after you open the foil pouch) on the actuator label using permanent ink.




2. Prime your AIRSUPRA inhaler before first use

Take the cover off the mouthpiece by squeezing both sides of the cover while pulling it down.



Shake the inhaler well. Hold your AIRSUPRA inhaler facing away from you and press down firmly on the top of the puff indicator to spray 1 test-puff into the air. Repeat the priming steps for a total of **4 test-puffs**. Shake your AIRSUPRA inhaler before each test-puff.

Shake before each test-puff. Repeat 4 times.



Extra puffs are provided for priming. **Do not skip priming.**

When and how to re-prime your AIRSUPRA inhaler


Re-prime your inhaler:

- if it is not used for more than 7 days.
- if it is dropped
- after rinsing white actuator

To re-prime your inhaler:

Take the cover off the mouthpiece and spray **2 test-puffs into the air**. Shake your AIRSUPRA inhaler before each test-puff.

Shake before each test-puff. Repeat 2 times.




Extra puffs are provided for priming. **Do not skip priming.**

How to use your inhaler – Inhaling your medicine

1. Check


Take the cover off the mouthpiece. **Check the mouthpiece for foreign objects and remove objects before use.**



Check inside mouthpiece.


1. Remove canister

Remove the canister and set it aside. Take the cover off the mouthpiece.




2. Shake and Inhale

Shake well




Breathe out fully



Place the mouthpiece into your mouth and close your lips around the mouthpiece.



Start to breathe in deeply and slowly while spraying **1 puff**. Continue breathing in until you cannot anymore.




Remove the mouthpiece from your mouth. Hold your breath for as long as you can or up to 10 seconds.




2. Rinse through both ends


Rinse each end of the actuator with warm water for 30 seconds. Shake off as much water as you can. **Do not use soap. Do not dry with a towel or tissue.**



60 seconds total



Look into the mouthpiece and actuator to make sure any medicine build-up has been completely washed away. If there is any build-up, repeat the weekly rinse step 2 again.



Check inside mouthpiece and actuator.

3. Repeat

Immediately repeat the shake and inhale step to take a second puff.



2 puffs equals 1 dose.


4. Place the mouthpiece cover back on

Place the mouthpiece cover back on the mouthpiece.




5. Rinse

Rinse your mouth with water, if available (**Do not swallow the water**).



3. Air-dry

Let the actuator air-dry, such as overnight.



Weekly rinse – Rinsing your actuator 1 time each week


- Rinse the white actuator 1-time each week so that the medicine does not build up and block the spray through the mouthpiece.
- **Do not allow the canister to get wet.**
- Re-prime your AIRSUPRA inhaler after rinsing.

4. Re-assemble and re-prime

Re-assemble AIRSUPRA when the actuator is fully dry.

First, place the mouthpiece cover back on. After the mouthpiece cover is on, gently press the canister down into actuator. **Do not insert the canister while the mouthpiece cover is off.**

Re-prime the inhaler by following the Re-prime steps above.



Emergency inhaler use when wet

If you need to use the inhaler before it is dry:

- shake off as much excess water as you can
- insert the canister into the actuator
- shake and spray 2-test puffs into the air away from your face
- take the prescribed dose
- rinse the actuator again following the Weekly Rinse steps.

Storing your inhaler

- Store your inhaler at room temperature between 68°F to 77°F (20°C to 25°C).
- **Do not store in a humid environment, such as a bathroom.**
- Keep your inhaler and all medicines out of the reach of children.

Throwing away your AIRSUPRA inhaler

Throw away your AIRSUPRA inhaler in the household trash:

- when the dose indicator reaches zero “0” or 6 months after you open the foil pouch, whichever comes first.

Do not reuse or use the actuator with medicine canisters from other inhalers.

Do not puncture or throw the canister into a fire or incinerator.

Ordering a new AIRSUPRA inhaler

- Order a new AIRSUPRA inhaler when the pointer on the dose indicator is in the yellow zone.

This Instructions for Use has been approved by the U.S. Food and Drug Administration.

Issued: January 2023

01/23 US-70802 1/23